



# CS 7200 Neo Edition IDS Information

**NEW**  
+ ETHERNET  
+ SMART PLATES  
+ SCAN & GO



# CS 7200 Neo Edition

## AT A GLANCE



Film-like imaging plates




19  
LP/MM

High resolution images



Most common sizes




Smarter Scan & Go workflow



8  
seconds

Fast image read-out





Powered with CS Imaging

# All common plate sizes

IDEAL FOR EVERYDAY NEEDS



- Adult periapical
- Bitewing



- Bitewing (children/small mouth)
- Can be used for periapical



- Periapical pediatric exams

# Speeding up your workflow

## FOUR EASY STEPS



**1** Identify the plate electronically



**2** Patient image is captured onto the plate.




**3** Plate is scanned and erased in one step. The plate is ready for reuse.



**4** Image appears within seconds, no click required.

# Scan & Go simplicity

## SMART WORKFLOW FOR GREATER EFFICIENCY



Plates are identified electronically prior to the exam

Plates are automatically recognized during scanning

Images are automatically sent to the appropriate computer and patient file

- Truly designed for multiple patients and multiple users
- Multiple users can use the system at the same time - without waiting in a queue
- Plates can be scanned in any order with no need for extra attention
- Eliminates risk of image mix ups and confusion

# The right accessories for the right results



x4

Smart Plates  
Size 0, 1 or 2 package  
4 plates per package



x200

Hygienic sheaths  
Size 0, 1 or 2 package  
200 pieces per package



Plate cleaning solution



Plate storage box and pad

# Which system is right for you


## QUICK COMPARISON



	CS 7200 Neo Edition	CS 7600
Use	Chairside / Shared use	Shared use / Multi-users
Theoretical resolution (max.)	25 lp/mm	25 lp/mm
True resolution (max.)	19 lp/mm	17 lp/mm
Smart Plates	✓	✓
Plate sizes	0-1-2	0-1-2-3-4
Examinations	Pediatric - Periapical - Bitewing	Pediatric - Periapical - Bitewing Long Bitewing - True occlusal
Time to first image	8 sec.	5 sec.
LCD screen with user instructions	-	✓

# Technical specifications

CS 7200	
Theoretical resolution (max)	25 lp/mm
Theoretical resolution (max)	19 lp/mm
Imaging plate sizes (Height x width)	Size 0 – 22 mm x 35 mm Size 1 – 24 mm x 40 mm Size 2 – 31 mm x 41 mm
Plate erasing	Automatic
Scan time	As little as 8 seconds for one image
Connectivity	USB and Ethernet
Power supply	100–240 V (ac), 50/60 Hz, 1.2 A
System dimensions (W x H x D)	130 x 270 x 300 mm 5.1 x 10.6 x 11.8 in.
Weight	3.5 kg





# Imaging plate systems portfolio

EXPLORE OUR HIGH-QUALITY IMAGING PLATE SYSTEMS

**NEW**  
**CS 7200 Neo Edition**  
*USB*  
Size 0,1,2  
**Your ideal chairside companion**

**NEW**  
**CS 7200 Neo Edition**  
*Ethernet Scan & Go*  
Size 0,1,2  
**Easy-to-use network reader for shared use**

**CS 7600**  
*Ethernet Scan & Go*  
Size 0,1,2,3,4  
**Versatile network reader for the widest range of diagnostic needs**

# CS 7200

## DELIVERY CONTENT

### Chairside configuration cat# 5942727

- CS 7200 scanner
- Insertion panel size 0, 1, 2
- Plates size #2 (4 pieces)
- Hygienic sheaths size #2 (200 pieces)
- Screen Cleaner sample bottle (plate cleaning solution)
- **USB cable**
- Power supply & cord
- International AC adapter
- Quick user guide (poster)
- User and installation guide (USB key)

### Network configuration cat# 5944608

- CS 7200 scanner
- **Scan & Go device (x2)**
- Insertion panel size 0, 1, 2
- Plates size #2 (4 pieces)
- Hygienic sheaths size #2 (200 pieces)
- Screen Cleaner sample bottle (plate cleaning solution)
- **Ethernet cable**
- Power supply & cord
- International AC adapter
- Quick user guide (poster)
- User and installation guide (USB key)



Carestream's Dental Imaging Software should be ordered separately  
Additional accessories can be ordered separately

# CAT Numbers

## EQUIPMENT

CAT #	Description	RRP
5942727	CS 7200 IMAGING PLATE SYSTEM	5999 € / £5499
<b>NEW CAT#</b> 5944608	CS 7200 + 2 x Scan & Go Devices	6499 € / £5999

## ACCESSORIES AND CONSUMABLES FOR CS 7200

CAT #	Description	RRP
5942511	CSD Smart Plate #0 - 4 Pcs	229 € / £197
5942537	CSD Smart Plate #1 - 4 Pcs	229 € / £197
5942545	CSD Smart Plate #2 - 4 Pcs	229 € / £197
5942503	Hygienic Sheaths #0 - 200 Pcs	39 € / £32
5942529	Hygienic Sheaths #1 - 200 Pcs	39 € / £32
5942578	Hygienic Sheaths #2 - 200 Pcs	39 € / £32
<b>NEW CAT#</b> 5944731	Scan & Go Device	279 € / £229
1029594	CS 7600 / CS 7200 Transportation Box and Pad	59 € / £54

# Timing and Preorders

- **Official Launch & deliveries start:** Mid April 2023
- **Launch event:** IDS, March 2023
- **Preorders:** accepted starting IDS with **special order forms** to avoid mixup between the current CS 7200 for which we still have stock and the CS 7200 Neo Edition, as well as current plates versus new plates.
  - [Export New Order Form - EN](#)
  - [Export New Order Form - FR](#)
  - [Export New Order Form - DE](#)
  - [Finished Goods Order Form - DE](#)
  - [Finished Goods Order Form - FR](#)
  - [Finished Goods Order Form - EN](#)
  - [Finished Goods Order Form - ES](#)
  - [Finished Goods Order Form - IT](#)