

Ergonomics, Access and Visibility

Ti-Max Z99L 1:5 achieves perfect Product Shopper score of 5.0



Scan/Click
For info article

BDIA DENTAL SHOWCASE
Raising the standards of oral care
13 - 14 March, 2026 • ExCeL London **STAND E50**



New NLZ E

Access meets endo

The unique motor that covers all aspects of endodontics – endodontics with a difference.



BUY **NLZ E** WITH **Ti-Max Z99L** FOR **£3,363**

SAVE £1666

See page 4 for more NLZ E offers →



New Ti-Max Z99L

New addition to the Ti-Max range of Speed Increasing Contra-Angles with a 100° angled head to enable easier access.



1 + 1 FREE* RRP £1,666*

See page 3 for more offers →

Restorative - Turbines

DENTAL SHOWCASE
 Raising the standards of oral care
 13 - 14 March, 2026 · ExCeL London
JOIN US AT STAND E50

Ti-Max Z HIGH PERFORMANCE TURBINES



4 years for iClave users



2 + 1 FREE* + 1 FREE Z micro



NEW

ULTIMATE POWER for quicker preps and cutting Zirconia

• Lightweight titanium body with Duragrip coating to reduce hand fatigue

Scan/Click
BUY NOW →



- **Optic**
- **ULTIMATE Power 30-44w****
- **Std head: Ø 12.1 x 12.7 mm**

MODEL Z990L
 CODE PA23730001
 RRP £1,214

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®, Bien-Air®

**At 3 - 4.2 bar



- **Optic**
- **Power 28-41w****
- **Mini head: Ø 10.6 x 10.8 mm**

MODEL Z890L
 CODE PA23720001
 RRP £1,214

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®, Bien-Air®



- **Optic**
- **Power 20w**
- **Micro head: Ø 9 x 10.8 mm**

MODEL Z micro (ultra-mini head)
 CODE PA23570001
 RRP £1,214

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®
 For regular or short shank burs

Scan/Click
 For full info on models/order codes



S-Max M



3 years for iClave users



OPTIC 3 + 1 FREE + 1 FREE Pico*



• Stainless steel • Enhanced power • 4 hole spray
 • DIY cartridge change • Clean head system

Scan/Click
BUY NOW →

NON-OPTIC 1 + 1 FREE*



- **Optic**
- **Power 26w**
- **Std head: Ø 12.1 x 13.3 mm**

MODEL M900L
 CODE P1254001
 RRP £800

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®, Bien-Air®



- **Optic**
- **Power 23w**
- **Mini head: Ø 10.6 x 12.4 mm**

MODEL M800L
 CODE P1253001
 RRP £800

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®, Bien-Air®



- **Optic**
- **Power 9w**
- **Ultra-mini head: Ø 8.6 x 9 mm**

MODEL Pico
 CODE P1140001
 RRP £1,160

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®, Bien-Air®

Used only with short shank burs.



- **Non-Optic**
- **Power 26w / 23w**
- **Std and mini heads**

MODEL M900
 CODE P1256001
 RRP £745

MODEL M800
 CODE P1255001
 RRP £745

M900 also available in fittings for KaVo®, W&H®

Pana-Max2



2 years for iClave users



1 + 1 FREE*

• Non-Optic turbine with single hole spray • Clean head system
 • Save time/money with DIY cartridge change • **12 month warranty**



- **Non-Optic**
- **Head: Ø 11.2 x 13.4 mm**
- **Coupling connection**

MODEL Pana-Max2 PTL
 CODE P1211001
 RRP £484

for NSK PTL coupling

MODEL Pana-Max2 KV
 CODE P1212001
 RRP £484

for KaVo® MULTIflex® coupling



- **Non-Optic**
- **Head: Ø 11.2 x 13.4 mm**
- **Direct connection**

MODEL Pana-Max2 M4
 CODE P1154001
 RRP £439

For Midwest M4 hose

MODEL Pana-Max2 B2
 CODE P1155001
 RRP £439

for Borden hose

Scan/Click
BUY NOW →

PANA SPRAY Plus

1 + 1 FREE*



Specially formulated handpiece and air motor lubricant.

MODEL PANA SPRAY Plus
 CODE Z182100
 RRP £35.20

Scan/Click
BUY NOW →

*Terms and conditions apply. Offers only available through participating dealers.

Restorative - Contra Angles

Ti-Max Z



4 years for iClave users



2 + 1 FREE + 1 FREE Z micro*



- Lightweight external titanium body • **Optic with internal water** • Small head/narrow neck
- **Unique spray switch on Z99L, Z95L, Z45L, Z25L**



- **Optic**
 - **1:5 Increasing**
 - Head size 08.9 x 13.4mm, at 100° angle for **improved access**
- MODEL Z99L
 CODE CA13140001
 RRP £1,666

NEW

- **Optic**
 - **1:5 Increasing**
 - Head size 08.9 x 13.4mm
- MODEL Z95L
 CODE C1034001-SW
 RRP £1,666

- **Optic**
 - **1:5 Increasing**
 - 10% shorter & lighter than Z95L
 - Head size 08.9 x 13.4mm
- MODEL nano95LS
 CODE C1099001
 RRP £1,586

- **Optic**
 - **1:4.2 Increasing + Spray switch**
 - 45° head
 - Head size 09.6 x 14.5mm
 - Ideal for minor oral surgery too
 - For regular or extra long (25mm) FG burs
- MODEL Z45L
 CODE C1064001
 RRP £1,666

- **Optic**
 - **1:1 Direct Drive**
- MODEL Z25L
 CODE C1038001-SW
 RRP £1,320

1 + 1 FREE*

S-Max M



3 years for iClave users



3 + 1 FREE + 1 FREE Pico*



- Stainless steel • **All with internal water** • Choice of optic & non-optic



- **Optic**
 - **1:4 Increasing**
 - Head size 9.9 x H.14.5mm
- MODEL M94L
 CODE CA1317001
 RRP £1,076

NEW

Delivers high torque for smooth and efficient removal of zirconia

- **Optic**
 - **1:5 Increasing**
 - Head size 9.9 x H.14.5mm
- MODEL M95L
 CODE C1023001
 RRP £1,076

- **Non-Optic**
 - **1:5 Increasing**
 - Head size 9.9 x H.14.5mm
- MODEL M95
 CODE C1026001
 RRP £893

- **Optic**
 - **1:1 Direct Drive**
- MODEL M25L
 CODE C1024001
 RRP £871

- **Non-Optic**
 - **1:1 Direct Drive**
- MODEL M25
 CODE C1027001
 RRP £572

- **Optic**
 - **4:1 Reducing**
- MODEL M15L
 CODE C1025001
 RRP £1,019

- **Non-Optic**
 - **4:1 Reducing**
- MODEL M15
 CODE C1028001
 RRP £681

- **Non-Optic**
 - **1:1 Straight**
- MODEL M65
 CODE H1008001
 RRP £559



FXm



2 years for iClave users



1 + 1 FREE*

- **Non-Optic Straight**
 - **1:1 Direct Drive**
 - **Push button chuck**
 - **No internal water**
- MODEL FX65m
 CODE H1014001
 RRP £255

- **Non-Optic Contra-angle**
 - **1:1 Direct Drive**
 - **Push button chuck**
 - **No internal water**
- MODEL FX25m
 CODE C1052001
 RRP £282

- **Non-Optic Prophy Contra-angle**
 - **4:1 speed reducing**
 - **Max speed 1,000 rpm**
 - **Brushes/Cups: Screw in**
- MODEL FX57m
 CODE C1054001
 RRP £317

FX

3 FOR £201, 10 FOR £568*

New modern version of the classic NAC-EC

- **Non-Optic Contra-angle**
 - **1:1 Direct Drive**
 - **Latch-type chuck**
 - **No internal water**
- MODEL FX22
 CODE C1129001
 RRP £165

NEW



*Terms and conditions apply. Offers only available through participating dealers.

Restorative - Electric Micromotors

NLZ E **NEW**

PORTABLE MICROMOTOR

Including reciprocating and rotary Endo features

NEW compact portable electric micromotor, for general & endodontic applications with contra angles (ie speed increasing/red band, 1:1 direct drive/blue band or speed reducing/green band).

- Lighter motor with enhanced power/torque
- Ideal for general (restorative, prophylaxis etc) and endo applications (rotary & reciprocating* modes)
- Enhanced patient safety function indicates when the handpiece needs repair to minimise it overheating
- Versatile installation options – portable or semi built-in

* The 'Reciprocating Endo Mode' can be operated only with the NLZ Endo Contra-angle and WAVE ONE® file. The NLZ requires access to a plug socket, and the use of an existing Midwest hose from the cart to gain access to coolant water/air and the existing foot control.

MODEL NLZ E
CODE Y1004004
RRP £3,363



BUY NLZ E WITH Ti-Max Z99L
FOR **£3,363** **SAVE £1666**



BUY NLZ E WITH Ti-Max Z99L
& Z25L FOR **£3,891** **SAVE £2458**



BUY NLZ E WITH
Ti-Max Z99L &
Z25L & NLZ ENDO
FOR **£4,037**
SAVE £2677



NLX

RETRO FITTED MICROMOTORS

Electric micromotors are the driving force behind speed increasing handpieces and can be fully integrated into the delivery system of most new and existing dental chairs. They can be retro fitted by an engineer or come as an option when buying new chairs. Controlling the specific speed and torque of micromotors within existing platforms is made easy thanks to built-in control pads.

Speed increasing contra-angles not only do everything a turbine does using the same FG burs, they provide a constant torque from the micromotor that does not decrease with resistance and higher load with no loss of power as they come in contact with the tooth. Higher load with no loss of power as they come in contact with the tooth.

To discuss the offers & options contact NSK:

Scan/Click



For more info



OFFER 1 BUY 3 x  AND GET **OPTION A** FREE



OFFER 2 BUY 4 x  AND GET **OPTION B** FREE



OFFER 3 BUY BUILT-IN OPTION A or B AND GET
FREE Ti-Max Z95L SPEED INCREASING CONTRA-ANGLE (1:5)



OPTION A NLX nano micromotor & BUILT-IN KIT*
(including MultiPad) RRP £2,387

OPTION B NLX plus micromotor & BUILT-IN KIT*
(includes Rotary Endo functions)
(including MultiPad) RRP £2,871

OPTION C Varios170 Lux ULTRASONIC HANDPIECE
(tip options for scaling, perio, endo & restorative procedures)
& BUILT-IN KIT* (including MultiPad)
RRP £1,553 (Y1004145)

OPTION D Option A + C
RRP £3,781

OPTION E Option B + C
RRP £4,238

BDIA DENTAL SHOWCASE

NSK ASPIRE Dental Academy WORKSHOPS

Workshops presented by **Raheel Malik** and **Richard Porter**



Anterior Veneer Preparation
Friday, 13 March | 11:30 - 13:00

Precision in Posterior Onlay Preparation
Saturday, 14 March | 10:00 - 11:30



Mobile dentistry

VIVA ace

MOBILE DENTISTRY AND SURGERY
BACK UP – ALL IN ONE SINGLE BOX

- Lightweight micromotor with additional endo functions
- Multipurpose ultrasonic handpiece
- Quick setup and storage
- 4-level adjustable height and angle



Width 462 mm, Height 284 mm, Depth 292 mm



MODEL VIVA ace
CODE Y1003773
RRP £13,408

Weight 8.6kg

FREE* 1 x



M95L



M25L



M65

Monthly RENTAL £433* Fully maintained for 3 years

Endodontics

NLZ E NEW

PORTABLE MICROMOTOR

Including reciprocating and rotary Endo features

NEW compact portable electric micromotor, ideal for endodontics.

- Lighter motor with enhanced power/torque
- Ideal for endodontic applications (rotary & reciprocating* modes)
- Enhanced patient safety function indicates when the handpiece needs repair to minimise overheating
- Versatile installation options – portable or semi built-in



MODEL NLZ E
CODE Y1004004
RRP £3,363

For OFFER see page 4

Reciprocating Endo Mode is operated with the NLZ Endo Contra-angle and WAVE ONE® file.

Varios Tips

ULTRASONIC TIPS FOR ENDO

Designed to aid agitation and irrigation. Compatible with stand alone Varios 970 / Varios Combi Pro, Varios 170 (built-in).



- E5 & E6 for lateral condensation
- E4D (diamond coated) for cleaning/smoothing or widening canals
- E7 / E8 for removing filling material from the canal
- E15D (diamond coated) for locating root canals & decalcification

MODEL Varios Endo tip kit
(tip stand + 6 tips above)
CODE Z221086/6

OFFER £264 Equivalent of 3+3 free tips + free tip stand

Additional tips RRP £85.90 Buy 1 get 1 free (diamond tips £92.20 Buy 1 get 1 free)



Scan/Click
For info article



Endo-Mate TC2

CORDLESS ENDO
MOTOR/HANDPIECE

The classic cordless endo motor/handpiece

- Use with rotary files

MODEL Endo-Mate TC2
CODE Y1001027
RRP £1,830



Endo-Mate DT2

CORDED ENDO
MOTOR/HANDPIECE

- Use with rotary files

MODEL Endo-Mate DT2
CODE Y1004230
RRP £1,730



iPex II

APEX
LOCATOR

MODEL iPex II
CODE Y1002208
RRP £1,086



BUY Endo-Mate DT2 AND GET iPex II Apex Locator FREE WORTH £1,086*

Powder Therapy

FLASH pearl

Supragingival powder for tough stains. Used with ProphyMate or the Varios Combi Pro (Prophy Chamber).

- Calcium carbonate powder (54µm)



4 x 300 g bottles
CODE Y900698
RRP £147.70

100 x 15 g sachets
CODE Y900693
RRP £201



SOFT pearl

Supragingival powder for moderate stains. Used with Prophy-Mate or the Varios Combi Pro (Prophy Chamber).

- Glycine based powder (65µm).
Water soluble, making clean up easier.



NEW

4 x 250 g bottles
CODE Y1500721
RRP £142.40



BUY ANY 2 PACKS OF FLASH OR SOFT pearl
GET A **FREE** nozzle for Varios Combi Pro or Prophy-Mate*

Perio Mate Powder

For subgingival biofilm removal or mild supragingival stains. Used with Perio-Mate or the Varios Combi Pro (Perio Chamber).

- Glycine based powder (25µm)



POWDER THERAPY PATIENT LITERATURE

Available on request with Prophy-Mate neo, Perio-Mate and Varios Combi Pro purchases.



BUY 2 PACKS OF Perio-Mate Powder
GET A **FREE** JET NOZZLE*


Ultrasonic Tips

Scan for full range of tips



1 + 1 FREE*

Scan/Click
BUY NOW



JIDA



JIDA SCIENCE



Read a clinical feature in JIDA by Pro Luigi Nibali and Dr Varkha Rattu including the use of NSK tips in Minimally Invasive Non-Surgical Therapy MINST

Just announced for June 2026

Prof Nibali and Dr Rattu will present **The Non-surgical & Surgical Periodontal Regeneration Course.**

Pre-course online theory session
Thursday, 4 June

Hands-on workshop
Saturday, 27 June



Scan for full details

Perio & oral hygiene

BDIA DENTAL SHOWCASE
Raising the standards of oral care
13 - 14 March, 2026 • ExCeL London
JOIN US AT STAND E50

Varios Combi Pro2 iPiezo engine®



NEW

AN EVOLUTION IN ULTRASONIC AND POWDER THERAPY

Gentle yet powerful, reliable and exceptionally durable, Varios Combi Pro2 has a host of advanced features designed to improve efficiency, control and enhance patient comfort. Engineered with meticulous attention to detail, Varios Combi Pro2 is an investment in clinical excellence and in the health of your patients.

- Superbly comfortable water temperature, for both ultrasonic & powder
- Durable titanium ultrasonic handpiece
- Jet Nozzle for Perio Handpiece
- Auto Cleaning
- Removable Powder Chambers*
- Multi-functional wireless foot control

Elevate your practice and patient care with the Varios Combi Pro2 - the future of ultrasonic and powder therapy.



MODEL Varios Combi Pro2
CODE Y1500505
RRP £6,750

RENT from £243pm†

Note: Connection to external 6mm air & water line sockets are needed to run this unit.

Call for a demo

Contact our team
Details on back cover

†Visit mysnk.co.uk/rental for more information.
*Comes with detachable Prophy Chamber Set.
Perio Chamber Set available as an optional extra.

FREE PERIO KIT2 & CLM MODULE SET WORTH OVER £1,400



Heated Water

Both ultrasonic and powder handpieces inject water that is heated to the approximate temperature of the oral cavity, ensuring greater comfort, even for those patients with sensitivity.



Jet Nozzle for Perio Handpiece

Executes precise, targeted treatment to areas with staining or thick biofilm on roots or restorations.

MODEL Jet nozzle for Perio Powder Handpiece
CODE TA10770106
RRP £107.60



Titanium Ultrasonic handpiece

- Durable & lightweight, with unrivaled ergonomics
- Non-slip & scratch-resistant Duragrip coating
- Slim LED handpiece enhances visibility

MODEL VA2-LUX-HP Ti
CODE EA1203X050
RRP £945



Watch the VCP2 demo with Cat Edney

iCart Prophy2

User-friendly and convenient trolley cart with the optional extra of a water bottle attachment, to create an independent water supply. This means that when the Varios Combi Pro2 is used with the iCart & bottle, you can treat patients without connecting to the external water mains supply.

MODEL iCart Prophy 2 (cart only)
CODE S9098
RRP £1,077

MODEL Water Bottle Kit
CODE Y1500719
RRP £725

MODEL Perio Kit2
CODE Y1500470
RRP £1,274



*Terms and conditions apply. Offers only available through participating dealers.

BDIA DENTAL SHOWCASE

NSK LECTURES



Dr Varkha Rattu
Minimally Invasive Therapies
Main Stage
Friday, 13 March



Cat Edney
Design and Use of Ultrasonic Tips
Oral Health Theatre
Saturday, 14 March

Perio & oral hygiene

Varios 970 iPiezo engine®

ULTRASONIC THERAPY

Ultrasonic scalers not only offer patients a more comfortable experience as less water is needed, but offer the user a wide range of 70 tips, not only for scaling/perio/implant maintenance, but for endo, extractions, restorative & MI applications.

Units come with:

- 1 handpiece
- G4, G8 & G16 tips
- 3 tip wrenches
- sterilisation case

• Optic

MODEL Varios 970 LUX
CODE Y1004161
RRP £3,103

• Non-Optic

MODEL Varios 970
CODE Y1004160
RRP £2,671

Scan/Click
BUY NOW



OFFER

FREE

2 x



2 x



Varios2 handpiece
• Optic

MODEL VA2-LUX-HP
CODE E351050
RRP £568



Varios2 handpiece
• Non-Optic

MODEL VA2-HP
CODE E350050
RRP £317

1 + 1 FREE*

Scan/Click
BUY NOW

1 + 1 FREE*

Scan/Click
BUY NOW

Air Scalers



Direct connection air scaler

- Oscillation frequency 5,600-5,850Hz
- Kit includes 3 scaler tips

MODEL AS2000
CODE T698001 (Midwest M4 hose)
CODE T697001 (Borden B2/B3 hose)
RRP £510 each



Optic air scaler with coupling connection

- Oscillation frequency 6,200-5,850Hz
- Scratch resistant with titanium body
- Kit includes 3 scaler tips

MODEL Ti-Max S970L
CODE T1005001
RRP £1,433

Also available in KaVo® fitting T1013001 and Sirona® fitting T1015001

FREE FX22
WORTH £165



FREE 2 x FX22 WORTH £330

TIPS



Scan/Click
Full range of tips

1 + 1 FREE*

Prophy Contra Angles



The design of the FX with its long slim neck (42mm vs 33.5mm) improves visibility of the operational field, even in difficult to access areas.



Scan/Click
BUY NOW



FX57m FX65m

- Max speed 1,000 rpm
 - Brushes/Cups: **Screw in**
- MODEL FX57m
CODE C1054001
RRP £317

- Max speed 40,000 rpm
 - Brushes/Cups: **Disposable Prophy heads**
- MODEL FX65m
CODE H1014001
RRP £255

1 + 1 FREE*



Scan/Click
BUY NOW

Scan/Click
BUY NOW

NEW

- Max speed 5,000 rpm
 - Brushes/Cups: **Snap on cups**
- MODEL AR-EC(k)
CODE Y110165
RRP £160

- Max speed 5,000 rpm
 - Brushes/Cups: **Screw in**
- MODEL AR-EC(s)
CODE Y110164
RRP £157

- **New modern version of the classic NAC-EC**
 - Max speed 30,000 rpm
 - Brushes/Cups: **Latex grip chuck**
- MODEL FX22
CODE C1129001
RRP £165

3 FOR £201, 10 FOR £568*

Perio & oral hygiene

BDIA DENTAL SHOWCASE
Raising the standards of oral care
13 - 14 March, 2026 • ExCeL London
JOIN US AT STAND E50

Prophy-Mate neo

SUPRAGINGIVAL POWDER THERAPY

- Compact, lightweight & portable
- Includes 60° and 80° nozzles
- 360° rotation
- Includes 6 sachets of FLASH pearl® powder

CODE Y135029 for NSK PTL coupling*
CODE Y135030 for KaVo® MULTiflex® LUX coupling
CODE Y135031 for Sirona® quick coupling
CODE Y135032 for W&H® Roto Quick® coupling
CODE Y135033 for Bien-Air® Unifix® coupling

RRP £1,203

Scan/Click
For more info



CHOICE OF OFFERS

1. FREE SECOND HANDPIECE WORTH £456* 

2. BUY WITH A Perio-Mate FOR £1,998* SAVE £576 

3. BUY WITH A iProphy CORDLESS PROPHY HANDPIECE FOR £1,425* SAVE £939 

PROPHY-MATE NEO NOZZLES

1 + 1 FREE*



CODE T797012 60° nozzle
RRP £128.50



CODE T797013 80° nozzle
RRP £128.50



IKIGAI

ORAL HYGIENE COMMUNITY

For all the latest in Perio & oral hygiene follow IKIGAI on

FACEBOOK

INSTAGRAM



Perio-Mate

SUBGINGIVAL POWDER THERAPY / BIOFILM REMOVAL

- Subgingival version of ProphyMate
- Eradicates biofilm in just 20 seconds
- Powder flow adjuster for delicate pockets
- Includes flexible plastic nozzle tips for deeper pockets
- Packed in a handy case with 160g powder
- Also used supragingivally on mild staining and around composite restorations/orthodontic brackets



Scan/Click
For more info



CODE Y1002654 for NSK coupling
CODE Y1002655 for KaVo® MULTiflex® LUX coupling
CODE Y1002656 for Sirona® quick coupling
CODE Y1002657 for W&H® Roto Quick® coupling
CODE Y1002658 for Bien-Air® Unifix® coupling

RRP £1,371

+ FREE SECOND HANDPIECE WORTH £456* 



Implantology & oral surgery

DENTAL SHOWCASE
 Raising the standards of oral care
 13 - 14 March, 2026 · ExCeL London
JOIN US AT STAND E50

Surgic Pro2

MICROMOTORS FOR IMPLANTS / ORAL SURGERY

Scan/Click
 For more about
 Surgic Pro2



Including the new stylish Surgic Pro2 with a cordless foot control, offering wireless connectivity to VarioSurg4 and Osseo 100+.

Complete Kits contain:

- Control unit
- Motor/Cord
- 20:1 Implant contra angle
- Foot control
- 2 year warranty

MODEL Surgic Pro2 (optic)
 CODE Y1004195
 RRP £6,429

MODEL Surgic Pro2 (non-optic)
 CODE Y1004196
 RRP £4,799



FREE



Osseo 100+

OR



iSD900

OR 2 x



20:1 implant
 contra-angles

ADDITIONAL MOTOR/CORD



• For Surgic Pro

MODEL SGL70M (optic)
 CODE E1023001
 RRP £2,854

MODEL SG70M (non-optic)
 CODE E1025001
 RRP £1,984

• For Surgic Pro2

MODEL SGL80M (optic)
 CODE Y1004211
 RRP £2,854

MODEL SG80M (non-optic)
 CODE Y1004212
 RRP £1,984

IRRIGATION TUBING



• Used with Surgic Pro / Pro2 / XT
 and VarioSurg 3 / 4



Irrigation tubing x5 (disposable)
 CODE Y900079
 RRP £66.50



Autoclavable Irrigation Tubing x1
 CODE S900055
 RRP £55.60



FREE 20:1 IMPLANT CONTRA-ANGLE OR SURGICAL STRAIGHT

WHEN YOU BUY AN ADDITIONAL MOTOR/CORD



SURGICAL STRAIGHTS 1:1

1 + 1 FREE*

- Used with standard HP burs
- **2 year warranty**
 (3yrs for iClave users)



• **Optic**
 MODEL Ti-Max X-SG65L
 CODE H1009001
 RRP £998

• **Non-Optic**
 MODEL Ti-Max X-SG65
 CODE H1038001
 RRP £681

IMPLANT CONTRA-ANGLES 20:1

1 + 1 FREE*

- Max speed: 2,000min⁻¹
- **2 year warranty** (3yrs for iClave users)



• **Optic**
 MODEL Ti-Max X-SG20L
 CODE C1003001
 RRP £1,331

• **Non-Optic**
 MODEL S-Max SG20
 CODE C1010001
 RRP £1,122



ZYGOMATIC CONTRA-ANGLE

1 + 1 FREE*

- **20:1 Reduction**
- **Twist chuck**
- **12 Month warranty** (2yrs for iClave users)



Handpiece for Zygomatic & Pterygoid implants. Its unique angulated shape offers optimal access and visibility.

• **Non-Optic**
 MODEL SGX-E20R
 CODE HA1200001
 RRP £1,637



*Terms and conditions apply. Offers only available through participating dealers.

VarioSurg 4 ULTRASONIC BONE SURGERY SYSTEM

NEW

FREE 12 EXTRA TIPS OR AN OSSEO 100+

With its stylish design and wireless foot control, VarioSurg 4 builds on the reliability and control of previous models to deliver a modern ultrasonic bone surgery system, with Piezo technology to minimise damage to surrounding tissue and bone.

Contains

- Slimline lightweight LED handpiece
- Stylishly designed control box
- Wireless foot control**
- Tip Stand with 6 tips - SG11, SG6D and SG7D for sinus lifts. H-SG1 for cutting bone with SG3 and SG5 for scraping.



MODEL VarioSurg 4
 CODES Y1500706 Full Kit (including cordless foot control)
 Y1500712 (without foot control)
 Y1500712 + Z1102003 (+ corded foot control)
 RRP All £8,661 each



NEW

SG12 Tip

NEW Ultra-slim tip 0.35mm, with blades on 3 sides for ridge splitting, chin bone harvests and general bone surgery, offering superior access from any angle.

CODE ZT02010001
 RRP £196



VarioSurg tips 1 + 1 FREE*



Scan/Click
 Extensive range of tips



Scan/Click
 For the core range of tips

LINK OPTIONS

Already have a Surgic Pro2?

Link it to a VarioSurg 4 and control both with 1 foot control**



VarioSurg 4 WITHOUT FOOT CONTROL - Y1500712

50% OFF FOR Surgic Pro2 USERS

+ FREE LinkStand 3 OR iCart Duo

Link Stand for VarioSurg 4 / Surgic Pro2

MODEL Link Stand3
 CODE ZA16230001
 RRP £321



iCart

MODEL iCart Duo
 CODE S9090
 RRP £1,077



Scan/Click
 Watch VarioSurg 4 and
 Surgic Pro2 link together



Osseo 100+

OSSEOINTEGRATION MONITORING DEVICE

Osseo 100+ measures implant stability and osseointegration to enhance decisions about when to load the implant. This is invaluable when using protocols with shorter treatment times and treating higher risk patients.

The simple steps to measure the ISQ value allow the implant loading period to be planned in advance.



30% OFF - ONLY £1,926*

MODEL Osseo 100+
 CODE Y1004176
 RRP £2,752



Scan/Click



Model: Multipeg
 RRP £44.70

More info &
 Multipeg List



“The Osseo device allows us to easily measure implant stability by assessing the increasing ISQ readings during the healing period. We use it routinely on sinus graft cases, when we place simultaneously, and also in cases with low stability at placement.”
The Campbell Clinic, Nottingham

*Terms and conditions apply. Offers only via participating dealers. **Only the Surgic Pro2 foot control (FC-81) is designed to link and control both units.

Decontamination & IPC

AFFORDABLE, RELIABLE & RAPID STERILISATION



iClave 18 & 24

B-type vacuum autoclave with 18 or 24 litre chambers



NEW

18 litre chamber
CODE DX91812B0100
RRP £6,637

24 litre chamber
CODE DX92412B0100
RRP £8,274

- Perfectly combines performance, user experience, & technology, ensuring an effective, comfortable, and worry-free experience.
- 5 trays of 2 different dimensions maximises capacity.
- High-efficiency steam generator allows for superior performance
- B-cycle 40min, S-cycle 31min
- User-friendly colour touchscreen
- Footprint: W451 x H492 x D669mm (smaller internal chamber on iClave 18)

iClave 18

FROM ONLY £3,995
+ install & delivery

RENT from £220_{pm}⁺
includes full maintenance package

iClave 24

FROM ONLY £4,995
+ install & delivery

RENT from £242_{pm}⁺
includes full maintenance package

iClave plus

B-type vacuum autoclave 18 litre



CODE DXY1003077
RRP £6,637

- Easy to use + USB data storage
- Fast 25min cycle for unwrapped loads
- 35min for pouched loads (+econ cycle 20min – one large tray; 15min – two small trays)
- 18L cylinder-shaped chamber offers 20% more space than traditional chambers & is ideal for implant kits
- Includes Helix Test + Strips
- Footprint: W445mm x H428 x D532mm

iClave plus

FROM ONLY £3,995
+ install & delivery

RENT for just £155_{pm}⁺

Sponsored by NSK



Find out more about the
**2026 IPC Team
of the Year Award**
on Stand G50 at
BDIA Dental Showcase.

12 MONTHS ADDITIONAL WARRANTY on new NSK handpieces or repairs for iClave users**

*Visit mynsk.co.uk/rental for more information. Every cycle over the agreed number of cycles at the end of the rental period incurs an additional cost of 7p per cycle - iClave 18 & 24 only.
**Additional warranty applies to NSK handpieces (new or repaired by NSK), that are still within their warranty period.

PROTECT IT

All-Inclusive Maintenance and Breakdown Contracts for NSK autoclaves from just **£64_{pm}**.

One annual maintenance visit or every 1500 cycles	Additional maintenance visit for units exceeding 1500 cycles	All service parts and service sundries included	All breakdown calls*	All breakdown parts including labour*	All engineer travel in the UK*	Free on-site operator training	Free technical software upgrades
---	--	---	----------------------	---------------------------------------	--------------------------------	--------------------------------	----------------------------------



More info on the NSK decon range and the Smart3 System of working

iCare

Handpiece lubrication unit



Precise dosing - less waste, better care

MODEL C2
CODE Y1002787
RRP £2,752

MODEL C3
CODE Y1002788
RRP £2,892

Connects to 6mm airline, 5-6 bar

NSK Eco Aware!



OFFER 1 BUY iCare C2 FOR £2,752 OR iCare C3 FOR £2,892 AND GET £2,137 WORTH OF INSTRUMENTS FREE

OFFER 2 BUY iCare FOR HALF PRICE C2 £1,376 OR C3 £1,446

RENT for just £73.95_{pm}⁺

CONTACT NSK

National Sales Manager
Mark Beckwith
mbeckwith@nsk-uk.com 07900 246529

Leics / Northants / CV & PE postcodes
Graham Brown
gbrown@nsk-uk.com 07541 864641

South East & East Anglia
Francis Faro
ffaro@nsk-uk.com 07900 246822

Scotland - Sales Manager
Angela Glasgow
aglasgow@nsk-uk.com 07525 911006

North West / North Wales
Clare Morris
cmorris@nsk-uk.com 07901 117669

South Wales
Robert Musgrove
rmusgrove@nsk-uk.com 07798 918741

Field Service Engineer - North
Jamie Duncombe
jduncombe@nsk-uk.com 07780 105457

N. Ireland & ROI
Asnain Sadiq
asadiq@nsk-uk.com 07435 840097

Derbys/Notts/Yorkshire/Lincs
Jo Tomlinson
jtomlinson@nsk-uk.com 07545 786619

South West
Ben Guest
bguest@nsk-uk.com 07350 361543

Field Service Engineer - South
Lee Morris
lmorris@nsk-uk.com 07435 783707

North East
Nicola Muirhead
nmuirhead@nsk-uk.com 07353 969563

West Midlands
Andy Hargreaves
ahargreaves@nsk-uk.com 07872 824175

South Central
Stephanie Collett
scollett@nsk-uk.com 07889 181099

NSK United Kingdom LTD www.mynsk.co.uk

Hertford House, Rutherford Close, Stevenage, Hertfordshire SG1 2EF, UK - info@nsk-uk.com - 0800 6341909



*Offers cannot be combined. NSK reserve the right to change prices and offers without notice. Published prices are RRP's and exclude VAT and installation. Offers only available through participating dealers. Additional warranties available from United Kingdom Ltd on new products and NSK United Kingdom Ltd repairs. To qualify for specific offers please send a copy invoice to NSK United Kingdom Ltd at Hertford House, Rutherford Close, Stevenage, Hertfordshire SG1 2EF, UK or email info@nsk-uk.com. Offers end 31/03/2026.