

#whdentalwerk



video.wh.com



Show casing W&H

Special offers
January - July 2026



Discover
Seethrough Flex
Augmented Reality

Turbines

High speed, low weight

All W&H turbines have a hygienic head¹ anti-retraction system to prevent cross infection.



Alegra TE-97/98

- › Without light
- › 3x Spray
- › Head Ø 10.4 / Ø 12.2 mm
- › 16 / 18 W
- › 390,000 / 330,000 rpm
- › Roto Quick® Connection

£ 265
RRP £586



Synea Fusion TG-97/98 - TG-97/98 L

- › Without light / LED+
- › 4x Spray
- › MINI Head Ø 10 / 11.5 mm
- › 18 W-390,000 rpm / 21 W-360,000 rpm
- › Roto Quick® Connection

TG-97/98 Non-optic
£ 519
RRP £1,037

TG-97/98 L Optic
£ 596
RRP £1,192

The TK-94 L microturbine is perfect for hard to access areas and the mouths of children.



Synea Vision TK-94 L

- › LED+
- › 5x Spray
- › MICRO Head Ø 9 mm
- › 14 W/410,000 rpm
- › Roto Quick® Connection

£ 774
RRP £1,548



Special Offer Price
When you buy
3x TG 97/98

x3



Synea Fusion TG-97/98 - TG-97/98 L

- › Without light / LED+
- › 4x Spray
- › MIDI Head Ø 10/11.5 mm
- › 18 W-390,000 rpm / 21 W-360,000 rpm
- › Roto Quick® Connection

3x TG-97/98 Non-optic
£ 1,399
RRP £3,111

3x TG-97/98 L Optic
£ 1,609
RRP £3,576

¹Exclusive W&H system patented in Europe and the USA



Air motors & Couplings

Flexible – for use with all standardised turbine tubing. Full compatibility with all turbines, air motors and scalers with Roto Quick and Multiflex connections.



360° instrument rotation for easy handling and to prevent the tube from twisting.



Couplings
Availables from
£ 112
RRP £217

Air motors
Availables from
£ 256
RRP £462

Contra-angles

High torque and precision



Synea Fusion WG-99 LT / WG-99 A

- › Contra-angle 1:5
- › Optic / Without light
- › 4x Spray
- › Head Ø 9.5 mm

£ 754
RRP £1,507

£ 674
RRP £1,348



Synea Fusion WG-56 LT / WG-56 A

- › Contra-angle 1:1
- › Optic / Without light
- › Single Spray
- › Head Ø 9.5 mm

£ 596
RRP £1,192

£ 476
RRP £951



Synea Fusion Power Edition WG-900 LT / WG-900 A

- › 1:4 high-speed contra-angle, designed to remove high-strength materials, like zirconia
- › Optic / Without light
- › 3x Spray
- › Head Ø 10.7 mm

£ 767
RRP £1,534

£ 688
RRP £1,375



Synea Vision Power Edition WK-900 LT

- › 1:4 high-speed contra-angle, designed to remove high-strength materials, like zirconia
- › Optic
- › 3x Spray
- › Head Ø 10.7 mm

£ 902
RRP £1,803



Alegra WE-56

- › Contra-angle 1:1
- › Without light
- › Optional External Spray
- › Head Ø 10.0 mm

£ 140
RRP £367



Alegra HE-43

- › Straight Handpiece 1:1
- › Without light
- › Optional External Spray

£ 130
RRP £319

Implantmed Classic SI-923



Implantmed Classic SI-923 Non-optic/Cable

- › Non-optic motor
- › Cable Foot Control
- › FREE WI-75 E/KM

£ 2,967

RRP £6,594

Oral Surgery & Implantology

Implantmed Plus SI-1023

Workflow driven solution



Implantmed Plus SI-1023 Kit 01 Non-optic/Cable

- › Non-optic motor
- › Cable Foot Control

› FREE WS-75

£ 3,421

RRP £7,602



Implantmed Plus SI-1023 Kit 02 Non-optic/Wireless

- › Non-optic motor
- › Wireless Foot Control

› FREE WS-75

£ 3,744

RRP £8,320

2 Years* Warranty W&H



W&H offer a synergistic solution for oral surgery guiding the clinician through the cycle from implementation to evaluation.



Implantmed Plus SI-1023 Kit 03 Optic/Cable

- › Optic motor
- › Cable Foot Control

› FREE WS-75 L

£ 3,951

RRP £8,781



Implantmed Plus SI-1023 Kit 04 Optic/Wireless

- › Optic motor
- › Wireless Foot Control

› FREE WS-75 L

£ 4,274

RRP £9,497

Free Surgical WS-75 Contra-angle Included

with selected Implantmeds

Scan function



Wireless foot control



Cable foot control



No light



Light



Complete Surgical Solution



Implantmed Plus SI-1023 + Piezomed Module Optic/Wireless

- › Optic motor and Wireless Foot Control
- › FREE WS-75 L
- › Piezomed Plus Module SA-435
- › Extraction Tip Set

£ 7,357

RRP £14,516

Implantmed Plus II

Smart assistance for advanced implantology



Implantmed Plus II Optic/Cable

- › Optic motor
- › Cable Foot Control

› FREE WS-75 L

£ 4,235

RRP £7,700



Implantmed Plus II Optic/Wireless

- › Optic motor
- › Wireless Foot Control

› FREE WS-75 L

£ 4,235

RRP £7,700



Implantmed Plus II Optic/Wireless/Scan

- › Optic motor
- › Wireless Foot Control
- › Scan function

› FREE WS-75 L

£ 4,648

RRP £8,450

New! In a Class of Its Own

Bundle Offer ImplantmedXpert Package

- › Implantmed Plus II SI-2102 Optic
- › Contra-angle WS-75 L
- › Piezomed Plus Module II SP-210

£ 7,739

RRP £15,279



Piezomed Classic Optic/Cable

- › Piezomed Classic
- › Handpiece with LED ring
- › Cable Foot Control
- › FREE BONE SET

£ 5,479

RRP £9,962



Piezomed Instruments



Full range of tips available



Bone surgery

£ 764

RRP £1,389



Extraction

£ 243

RRP £441



Lateral sinus lift

£ 549

RRP £997



Piezomed Plus Module

- › Power 24W. Handpiece with light SA-40 L
- › Automatic instrument detection
- › Instrument changer and instrument library
- › FREE EXTRACTION SET INCLUDED

£ 3,513

RRP £6,829



For Implantmed SI-1023 only



Piezomed Plus Module II

- › Power 24W. Handpiece with light SA-40 L
- › Automatic instrument detection
- › Instrument changer and instrument library
- › FREE EXTRACTION SET INCLUDED

£ 3,513

RRP £6,829



For Implantmed Plus II only



New!

Surgical Handpieces

High-quality stainless steel with a **scratch-resistant coating** makes W&H surgical handpieces robust and gives them a particularly long service life.



Self-generating LED Surgical Handpieces available.



Contra-angle 20:1 WI-75 E/KM

- › Without light
- › External single spray
- › Not dismantlable

£ 670

RRP £1,339



Contra-angle 20:1 WS-75

- › Without light
- › External single spray
- › Dismantlable

£ 828

RRP £1,655



Contra-angle 20:1 WS-75 L

- › Mini LED+
- › External single spray
- › Dismantlable

£ 905

RRP £1,810



Contra-angle WS-91 / 1:2.7

- › Without light
- › 45° head with 3x spray
- › Dismantlable

£ 780

RRP £1,559



Contra-angle WS-91 L / 1:2.7

- › Mini LED+
- › 45° head with 3x spray
- › Dismantlable

£ 863

RRP £1,725



Straight handpiece S-11 / 1:1

- › Without light
- › External single spray
- › Dismantlable

£ 688

RRP £1,336



Straight handpiece S-11 L / 1:1

- › Mini LED+
- › External single spray
- › Dismantlable

£ 941

RRP £1,881



Handpiece 20:1 Zygoma SZ-75

- › Without light
- › External single spray
- › Dismantlable

£ 1,341

RRP £2,682



Osstell Beacon

Measurement of implant stability



Osstell SmartPegs

The SmartPeg measures the resonance frequency when it is attached to an implant.

£ 88

RRP £121

£ 1,700

RRP £5,436

When you buy an Osstell Beacon register online at Osstell Connect to claim your 5 free* Smartpegs

*Free Smartpegs only available with the purchase of an Osstell Beacon and on registering for Osstell Connect

Osstell Classic

Osstell Classic helps you to objectively determine implant stability and measure the degree of osseointegration.


£ 1,390

RRP £2,884



Osstell Beacon

Measure and evaluate the stability of the implant and control its osseointegration quickly and easily thanks to the easy to read colour light indicator and ISQ reading. Reduce treatment times and avoid rejection in high-risk patients.



What are the positive outcomes when using the Osstell Implant stability solutions?

- › A positive experience for patients, encouraging trust through understanding
- › Enhanced patient education builds confidence, resulting in increased treatment acceptance
- › Reduction in treatment visits and failures with improved average healing times
- › Improved outcomes for patients with risk factors
- › Enhanced communication between specialists ensuring a more synergistic approach to treatment
- › Improved treatment documenting

Osstell ISQ devices with the Osstell ISQ Scale help implantologists to make informed implant treatment decisions.



Osstell Connect is designed to document quantified data which ensures traceability, whilst complementing data in your Patient Management System.



Dr Aly Virani speaks about how W&H streamlines his workflows:

Since I began using the W&H Piezomed unit, I've noticed that surgical procedures that involve bone removal, block bone grafting and sinus surgery have become safer and more predictable. Patients also report less postoperative pain versus traditional methods. Using the Piezomed alongside the Implantmed give me the flexibility I need during treatment. I made the decision to invest in Osstell Beacon because of the huge evidence base behind the technology. I knew that we needed a better way to monitor implants and to know when they are ready to be restored, what I didn't realise was the effect that this would have on our patients and our ability to communicate with them and bring them along for the journey. What we've seen is a huge amount of engagement and understanding from our patients to their body's response to implants and the occasions where more time is needed. What it has meant is that they see the ISQ rating, they see the scale, they see the technology and the evidence base and they really understand that we are trying to help them and that we are leading with the latest advancements. What this has meant for the practice is that we are sure on how to proceed and that we are able to communicate this to our patients, so I couldn't recommend the Osstell Beacon more highly.

The W&H Surgical Solution has been a real practice builder and led to fewer complications which is so important.

Dr Aly Virani North Cardiff Dental & Implants

Cleaning & Disinfection

BePro Professional Cleaning and Disinfection

Soaking bath

Soaking bath for the immersion of dental instruments in the decontamination phase.

£ 40
RRP £48



Disinfectant I

Powerful disinfectant concentrate for medical and dental instruments.

£ 41
Bottle 2l RRP £50



Disinfectant R

Ready-to-use solution for the disinfection and cleaning of rotary dental instruments.

£ 14
Bottle 1l RRP £17



£ 45
Refill 5l RRP £55

Disinfectant Wipes

120 wipes for the disinfection of non-invasive medical devices.

£ 10
Tube RRP £12

Disinfectant Spray

Surface disinfectant for the rapid disinfection and cleaning of non-invasive medical devices.

£ 14
Spray 1l RRP £17

£ 48
Refill 5l RRP £58

Sterilization, Hygiene & Maintenance

Internal cleaning & lubrication

Increase the life of your handpieces and reduce the costs of repairs!



10" quick connect by W&H

Assistina Twin

- Two process chambers allowing loading of a 2nd handpiece while the 1st is processed
- Superior Maintenance & Lubrication
- Intuitive and easy to use
- A flexible range of adaptors available

£ 1,484
RRP £3,299

Handpiece Maintenance Solution



Special Offer
**Lisa Mini 5 liters &
Assistina One**

£ 4,420
RRP £7,449



18" quick connect by W&H

Assistina One

Internal cleaning & lubrication of one handpiece in only 18 seconds, with Assistina One the aerosols stay in the HEPA filter thereby releasing clean air.

£ 1,050
RRP £2,235



U-son

- Available with 3 or 6 litre tank
- Stainless steel chamber
- Programmable temperature and timer
- Automatic 'de-gas' function

£ 535
3 liters RRP £824

£ 837
6 liters RRP £1,288

Pro Service

PremiumCare

Gives you the flexibility to select a plan that best suits the individual needs of your practice.

The standard PremiumCare package offers a flexible plan allowing you to add optional extras as you please or the **PremiumCare Plus or Flex**** plan, which include fully comprehensive service packages and the required pressure vessel inspection and validation testing at no extra cost.

PremiumCare

Maintenance plan for all W&H sterilizer models, with a maximum limit of 4,000 cycles over the 3 year term*.

PremiumCare Plus & Flex**

A fully comprehensive maintenance plan for all W&H sterilizer models, including PVI certification, unlimited cycles and validation in line with manufacturer's or current HTM 01-05 guidance (model dependent).

Our Technical Service Team
are here to support you nationwide!



W&H
ProService

Services	Premium Care	Premium Care Plus & Flex**
All routine service visits for three years	✓	✓
All service parts and scheduled sundries	✓	✓
All breakdown calls	✓	✓
All breakdown parts including labour	✓	✓
All engineer travel in the UK Excluding Highlands, Islands and ferries etc,	✓	✓
Free on-site operator training	✓	✓
Free technical software upgrades	✓	✓
Unlimited Cycles - for the duration of the plan	⊖	✓
PVI Certification - in line with Pressure Systems Safety Regulations 2000:128	⊖	✓
Periodic Test Logbook	⊖	✓
Validation in line with manufacturer's or current HTM 01-05 guidance (model dependent)	⊗	✓
Installation and Commissioning in line with current guidelines	⊗	✓
Helix Test Kit***	⊗	✓
Extra Discount	⊗	✓

✓ Included ⊗ Not Included ⊖ Optional

*An additional cycle tariff will be levied if the 4,000 cycle limit is exceeded within the 3 year period ** PremiumCare Flex Plans are only available for certain models *** Only applies to B Type Sterilizers



Thermal washer disinfector

The only system with two washing levels, for powerful cleaning, full flexibility, and faster turnaround in even the busiest clinics.



W&H
ProService
Pioneering
Excellence

Thed

- › Two washing levels
- › Fast cycles with high cleaning performance
- › Active air-drying with HEPA filter
- › Full connectivity options
- › Internal positioning of process chemicals
- › Unique cartridge system (tabletop)
- › Modern design
- › Comprehensive range of accessories
- › User-friendly interface

Offers include upper and lower rack

Thed Tabletop

£ 5,887
RRP £9,128

Thed Underbench

£ 6,505
RRP £9,746

Thed+

- Includes Thed features, plus:
- › EliSense: device status indicator
 - › Automatic rack recognition
 - › Automatic manifold injector recognition
 - › Automatic program recommendation
 - › Automatic adjustment of cycle time, water, process chemicals, and power consumption
 - › Washing chamber with internal light

Offers include upper and lower rack

Thed+ Tabletop

£ 6,887
RRP £9,850

Thed+ Underbench

£ 7,505
RRP £10,468



Cleaner TE

Tri-enzymatic detergent with neutral pH. No neutraliser required!

£ 47

Canister 5l RRP £58



Rinse Aid

For improved drying phase.

£ 19

Bottle 1l RRP £23

Sterilization at its Best!

With a wide range of sterilizers, W&H is able to offer the perfect option to meet individual needs.



Lyla S type Sterilizer New concept sterilizer

- › Mobile App available via upgrade for remote data storage
- › Future-proof upgradability to a type B sterilizer using the new Activation Code concept

- Unique upgrade flexibility with **Activation Codes**
- › B Cycle
- › Remote data storage
- › Traceability
- › All in One

£ 3,504

17 l. RRP £4,905

£ 3,915

22 l. RRP £5,159



Lara B type Sterilizer

- › Mobile App for Real Time Remote Monitoring
- › Eco Dry technology automatically adapts the drying time to the mass of the load
- › Unique optional upgradability to increase: Performance and Traceability or to add a Fast Cycle

- Unique upgrade flexibility with **Activation Codes**
- › Performance**
- › Traceability
- › Fast Cycle
- › All in One
- ** Only for Lara 17 and 22 litres

£ 4,790

17 l. RRP 7,063

£ 5,408

22 l. RRP £7,931

£ 5,923

Lara XL 28 l. RRP £8,747

£ 7,417

Lara XXL 38 l. RRP £10,682

Available in
17-38 Litre
versions

All W&H Sterilizers offer:

- › Intuitive user interface and touchscreen display for a smoother user experience
- › Cycle traceability available thanks to Wi-fi connection and built in data logger
- › Ergonomic design, with modular feet for flexible positioning
- › Easy to clean surfaces and accessible water
- › Automatic Water Fill function
- › Pre-programmable vacuum and helix (for B Type Sterilizers) test function to run automatically out of surgery hours



Lisa B type sterilizer

- › Full extended traceability thanks to EliTrace
- › Mobile App for Real Time Remote Monitoring
- › Eco Dry+ technology automatically adapts the drying time to the mass of the load
- › Elisense user interface generates status, temperature and other cycle information
- › Cycle times from 13 minutes for specific loads
- › State of the Art Vacuum B Sterilizer Technology for efficient, reliable, time saving, intuitive and ease of use

£ 5,717

17 l. RRP £8,722

£ 6,479

22 l. RRP £9,868



NEW Lisa MINI B type sterilizer

- › Cycle times from 10 minutes
- › Integrated vacuum pump
- › Small and easily stackable: W 44cm x H 20cm x D 67cm
- › Built in water tank and automatic water filling - no external storage
- › For up to 22 instruments, offers complete traceability
- › 5 l. model includes internal vacuum pump
- › Includes integrated water tank

£ 3,502

5 l. RRP £5,214

3 liters version available
Requires compressed air

Lisa Safe Label Printer & QR Barcode Reader

Quick and easy tool for linking the sterilization cycle and dental products used for treatment

Barcode reader

£ 150

RRP £242

Lisa Safe printer

£ 538

RRP £827



*with fast cycle upgrade

Sales Team

Get in touch with us!
We will treat your request immediately.



Tom James

National Sales and
Marketing Manager

07876 527720

tom.james@wh.com



Richard Mansbridge

Area: London, South East, Anglia,
East Midlands and Humberside

07775 792628

richard.mansbridge@wh.com



David Barrass

Area: Scotland, Humberside,
North of England and North West

07795 315089

david.barrass@wh.com



Caroline Caine

Area: West Midlands, Manchester, Merseyside
South West England, Isle of Man, Channel
Islands, Wales, Northern Ireland and Eire

07786 515331

caroline.caine@wh.com

To place an order contact W&H UK on

+44 (0) 1727 874990

or contact your preferred dental dealer

W&H (UK) Ltd

6 Stroud Wood Business Centre,
Park Street, St Albans, AL2 2NJ

t +44 (0) 1727 874990

f +44 (0) 1727 874628

office.uk@wh.com

www.wh.com



New products

Exciting news across the full
range of W&H products.

seethrough



implantmed^{PLUS}

Offers valid from: 1st January to 30th June 2026. Information is correct at time of publication and relates to W&H products and services available in the UK and Ireland. Prices are quoted in £ sterling and exclude VAT. Savings are based on RRP and rounded to the nearest £. Free goods are not transferrable and no cash equivalent is available. Terms and Conditions apply. E. & O. E.